Histoacryl®/Histoacryl® L

BG	Инструкции за употреба	– 2	LV	Lietošanas instrukcija	– 101
CS	Návod k použití	– 11	NL	Gebruiksaanwijzing	– 108
DA	Anvendelsesinstruktioner	– 18	NO	Bruksanvisning	– 116
DE	Gebrauchsanweisung	– 25	PL	Instrukcja użytkowania	– 123
EL	Οδηγίες χρήσης	– 33	PT	Informação para o utilizador	– 131
EN	Instructions for use	– 42	RO	Instrucțiuni de utilizare	– 139
ES	Instrucciones de uso	– 49	RU	Инструкция по применению	– 146
ET	Kasutusjuhend	– 56	SK	Návod na použitie	– 155
FI	Käyttöohje	– 63	SL	Navodila za uporabo	– 162
FR	Consignes d'utilisation	– 71	SV	Användarinstruktioner	– 169
HU	Használati útmutató	– 79	TR	Kullanım talimatı	– 176
IT	Istruzioni per l'uso	– 86	UK	Інструкція до застосування	– 183
LT	Naudojimo instrukcijos	– 94			

CE 0123

eIFU v.01 506675 - IFU(506530 02-04/23)

B | BRAUN

SHARING EXPERTISE

B. Braun Surgical, S.A.
Carretera de Terrassa, 121
08191 Rubí, Spain

Phone +34 93 586 62 00
Fax +34 93 586 10 96

www.bbraun.com

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EN

Description of the device

Histoacryl® is a tissue adhesive whose intended purposes are holding approximated tissue for closure of skin wounds, embolic agent for sclerotherapy of large oesophageal or fundal varices and fixating agent for fixation of hernia meshes.

Histoacryl® is a sterile, liquid tissue adhesive consisting of n-butyl-2-cyanoacrylate. For easy visual assessment of the thickness of the layer which has been applied, Histoacryl® is coloured with D&C violet No. 2 dye. A translucent version (Histoacryl® L), which does not contain dye, is also available.

Histoacryl® rapidly polymerises into a solid when in contact with body fluids, living tissues or blood.

Such Histoacryl® polymerised layer gradually and naturally detaches after some days when used for skin closure, or slowly degrades in the body when used for sclerotherapy and mesh fixation.

The clinical benefits of using Histoacryl® include wound closure and microbial barrier in skin closure, achievement of final haemostasis in sclerotherapy, and atraumatic fixation in mesh fixation.

Histoacryl® is supplied in 0.5 mL single use plastic ampoules.

Each ampoule is packed in an aluminium pouch to protect the ampoule from external contamination, keeping at the same time the external sterility of the ampoule.

Indications for use

1. **SKIN CLOSURE:** Closure of skin wounds without tension (including clean surgical incisions and incisions from minimally invasive surgery), and simple, thoroughly cleansed, trauma-induced lacerations.

2. **SCLEROTHERAPY:** Sclerotherapy of large oesophageal or fundal varices.

3. **MESH FIXATION:** Fixation of hernia meshes, especially in inguinal hernia surgery

Histoacryl® is intended to be used in different patient populations as follows:

- In adults, Histoacryl® is indicated for skin closure, sclerotherapy of large oesophageal or fundal varices and fixation of hernia meshes, especially in

inguinal hernia surgery.

- In children and neonates, Histoacryl® is only indicated for skin closure.

- Intended users are healthcare professionals that are qualified and familiar with the technique and with the properties of liquid tissue adhesives, specifically cyanoacrylates, before applying Histoacryl®.

Contraindications

- Histoacryl® is not indicated for holding together wound surfaces on: internal organs, the brain surface, central or peripheral nervous system, as such action could cause tissue damage and scarring, with the resultant complications.

- Histoacryl® is not to be applied to mucosal surfaces, mucocutaneous junctions, areas with dense natural hair or within the conjunctival sac of the eye.

- Histoacryl® is not to be applied on the surface of the eye. If Histoacryl® accidentally drops in the eye, adhesion can be avoided by immediately rinsing with water. If adhesion or substance depositions occur, especially at the eye, they will dissolve without further action in few days. Do not force the eye opening.

- Histoacryl® is not to be applied to dirty wounds or to wounds showing evidence of infection or gangrene, or to wounds of decubitus aetiology.

- Histoacryl® is not to be used on patients with known preoperative systemic infections, uncontrolled diabetes, or diseases or conditions that are known to interfere with the wound healing process.

- Histoacryl® has not been evaluated in patients with a history of hypertrophic scarring or keloid formation.

- Histoacryl® is not indicated on patients with a known hypersensitivity to cyanoacrylate, formaldehyde or the dye D&C Violet No. 2.

- In sclerotherapy, Histoacryl® is not indicated for the use in patients with congenital abnormalities leading to an increased risk of systemic embolisation.

- In sclerotherapy, Histoacryl® is not to be applied in those varices that have not ever bled.

- No more than 1.5 mL per treatment shall be applied in adult patients and no more than 0.5 mL per treatment shall be applied in paediatric patients weighing less than 10 kg.

No more than a total of 5 mL shall be applied to a single patient during her/his lifetime.

- Histoacryl® is not indicated for mesh fixation in IPOM procedures, given the risk of viscera adhesion.
- Histoacryl® is not indicated in paediatric population in sclerotherapy and mesh fixation.
- Histoacryl® is not indicated in pregnant, breastfeeding or women planning on pregnancy.

1. SKIN CLOSURE

Mode of action

Histoacryl® polymerises into a solid substance that attaches to the tissue in contact with tissue. Polymerisation is triggered by weak bases and nucleophilic functional groups commonly found in body fluids, living tissues or blood.

Studies(*) have shown that Histoacryl® acts as a barrier to microbial penetration as long as the adhesive film remains intact.

(*)Histoacryl® microbial barrier was in vitro tested during 7 days for the following bacteria: *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Escherichia coli*, *Pseudomonas aeruginosa*, *Enterococcus faecium*, *Brevundimonas diminuta* and *Candida albicans*.

Mode of application

- Prior to the application of Histoacryl®, provide local anaesthesia, if necessary, ensure that the edges of the wound have been thoroughly cleaned and debrided and always assure proper haemostasis. The areas to be joined shall be as dry as possible. Ensure that any remaining topical medications/anaesthetics are removed prior to application.
- Histoacryl® must be used in combination with and not in substitution of subcuticular sutures. Appose skin edges avoiding irregularities on the wound surface.
- Pull open the aluminium pouch and remove the ampoule containing the adhesive in sterile condition. The ampoule containing the adhesive shall only be removed from the aluminium pouch immediately prior to use.
- Hold the ampoule with the cannula pointing upwards and flick it sharply to ensure no adhesive remains trapped within the cannula. Open the ampoule

by twisting off the ribbed tip of the cannula keeping it until the end of the procedure. Hold the yellow part of the ampoule between two fingers, keeping it vertical with the thin end pointing upwards while opening the ampoule. This will prevent any loss of Histoacryl® by spillage from the ampoule when opening.

- Apply Histoacryl® to the skin surface by pressing the side walls of the plastic ampoule. Hold the edges together for approximately 30 seconds after application to allow Histoacryl® to cure and to prevent displacement of the wound edges.

- Avoid heavy application. Unless otherwise prescribed, as less possible Histoacryl® skin adhesive should be applied; the amount applied is sufficient if slight colouring is visible. A sparing application of thin layers is required for undisturbed wound healing.

- Discard any remaining Histoacryl® without contaminating the wound surroundings. Prior to discarding the ampoule, put back the ribbed tip that was removed when opening the ampoule.

- If removal of the polymerised Histoacryl® layer is necessary for any reason, carefully apply petroleum jelly or acetone to the Histoacryl® film to help loosen the bond. Peel off the film, do not pull the skin apart.

Warnings

- Histoacryl® must not be introduced into the wound since this would interfere with wound healing. This could also lead to long-term blue staining of the tissue. For facial application it is recommended to use Histoacryl® translucent.

- A heavy application may cause thermal damage to tissues, resulting in delayed healing. Connective tissue healing can be impeded when too much tissue adhesive is applied.

- Histoacryl® shall not be applied to wet or bleeding wounds. Excess of hydroxyl ions or nucleophilic groups usually found in alcohols, body fluids, living tissue and blood may accelerate polymerisation, resulting in the release of an excess of polymerisation heat and release of smoke.

- Histoacryl® is not to be applied subcutaneously since it might result in transient foreign body reaction that may occasionally take the form of an inflammatory

reaction.

- Histoacryl® is not intended to be applied to wounds subject to high skin tension or to areas of increased skin tension such as the elbows, knees, or knuckles.
- Histoacryl® is not to be used in areas of skin excision.
- Do not reuse. Reuse might cause an infection hazard for patients and/or users, impaired product functionality, and risk of injury, illness or death due to contamination.
- Histoacryl® is made from n-butyl-2-cyanoacrylate (NBCA) and D&C violet No. 2 dye in the coloured version, which are non-metallic and non-conducting materials, therefore it is considered MR safe.

Precautions

- Only product that is liquid in consistency must be used. The condition of the adhesive can be judged with the plastic ampoule unopened.
- The areas to be joined shall be as dry as possible. Excess hydroxyl ions or nucleophilic groups usually found in alcohols, body fluids, living tissue and blood may accelerate Histoacryl® polymerisation and might affect wound closure results. Any remaining topical medications/anaesthetics must be removed prior to application.
- If too much adhesive is accidentally applied, it can be removed during the first few seconds using a dry swab. Care must be taken; the dry swab could adhere.
- The skin surfaces should be held together for approximately 30 seconds after application of Histoacryl®.
- Wounds shall be kept dry following closure with Histoacryl®. Do not apply liquid or topical medications onto wounds closed with Histoacryl®. The polymerised film might be weakened, leading to skin edge separation.
- Histoacryl® should not be used on wound sites that will be subjected to repeated or prolonged moisture, movement or friction.
- Accidental bonding of unwanted skin may occur. Do not pull apart skin. If necessary, adhesions in less sensitive skin areas can be carefully dissolved with acetone or by soaking in warm water until the skin can

be separated.

- If using on the eye area, adjust the patient position to avoid accidental flow of Histoacryl® to unintended areas of the body. Use petroleum jelly to create a barrier close to the wound and use surgical gauze to protect the patient's eye.
- Care must be taken to avoid direct contact of instruments, cloths or gloves with the adhesive applied, as they might adhere to the tissue. Instruments stained with Histoacryl® can be cleaned with dimethyl formamide or acetone.
- As general practice, contact of the cannula with the skin should be minimized.
- Apply Histoacryl® as briefly as possible to reduce the risk of accidental bonding of the cannula to the skin.

Side effects

The use of this product leads to a mild exothermic polymerisation reaction. The released heat may result in a sensation of warmth. Histoacryl® migration below the epidermal surface may result in "tattooing" of the underlying tissue.

The following adverse events have been reported for Histoacryl®: wound dehiscence, wound infection, surgical site infection, bleeding, contact dermatitis, rash, bruising, blister formation and pruritus.

2. SCLEROTHERAPY

Mode of action

When instilled into the varix, Histoacryl® polymerises and hardens to a solid material, thereby, it seals the opening of the varix and occludes the blood flow into the varix, causing immediate cessation of bleeding (if active bleeding varix).

Within the varix, and following injection of Histoacryl®, there is an inflammatory reaction with foreign body thrombus formation, lymphohistiocytic infiltration and some foci of vascular wall necrosis. Finally there is a fibrosis in the region of the injection.

Mode of application

a) Pull open the aluminium pouch and remove the ampoule containing the adhesive in sterile condition. The ampoule containing the adhesive should only be removed from the aluminium pouch immediately prior to use.

b) Hold the Histoacryl® ampoule with the cannula pointing upwards and flick it sharply to ensure no adhesive remains trapped within the cannula. Open the ampoule by twisting off the ribbed tip of the cannula keeping it until the end of the procedure. Hold the yellow part of the ampoule between two fingers, keeping it vertical with the thin end pointing upwards while opening the ampoule. This will prevent any loss of Histoacryl® by spillage from the ampoule when opening.

c) Sclerotherapy with Histoacryl® in oesophageal or fundal varices, shall be only performed by a physician who is familiar with its application for this indication. The method of application of Histoacryl® depends on the user's individual experience in endoscopic varix therapy with cyanoacrylate adhesives.

d) Histoacryl® could be diluted by the addition of Ethiodized oil in order to adjust the polymerisation behaviour, with a proportion based on medical criteria. Values of *in vitro* polymerization of Histoacryl®/Ethiodized oil mixtures in blood can be found in Table 1 (caution: values obtained in static conditions, for information only. The conditions *in vitro* are not rigorously predictive of clinical conditions).

Table 1: Mixture proportions

Ratio (Histoacryl®: Ethiodized oil)	Polymerization time (s)
1:1	11.0 ± 3.8
1:5	39.3 ± 3.4

e) Immediately before and after the injection of the Histoacryl® mixture, the injection needle needs to be flushed with 5 % non-ionic glucose or dextrose solution to completely remove any anionic remnants from its lumen. This helps to prevent unintended polymerization of the mixtures by premature contact with anions.

f) To prevent premature solidification at the tip of the needle, it should not be primed with glue before injection but it might be soaked with Ethiodized oil.

g) Puncture the varix with the injection needle and

inject 1 mL maximum of the mixture intravariceal.

To minimize the risk of embolism, do not inject more than 1 mL of the mixture into the varix.

h) Complete obliteration is defined as an absolute firmness of the injected varix. Typically, if obliteration cannot be reached by one injection, additional injections at other site of the varix with 1 mL of the mixture should be used to stop the bleeding. The total amount of injections will be dictated by medical criteria considering the total dose limit specified in the *Contraindications* section.

Warnings

- The application of Histoacryl® must be limited only to large oesophageal or fundal varices. In case of small-calibre varices, administration must be carried intravariceal due to the fact that application of Histoacryl® outside the varix (paravariceal injection), can lead to necrosis of the stomach wall, deep ulceration, and perforation.

- The use of Histoacryl® induces an eradication of varices. Nevertheless, for some weeks there is a risk of recurrent bleeding, until the varices are completely thrombosed. Repeated Histoacryl® applications can be performed following medical criteria. The total dose limit specified in the *Contraindications* section must be considered.

- Patients with advanced liver disease, including ascites, have a higher risk of re-bleeding.

- Do not reuse. Reuse might cause an infection hazard for patients and/or users, impaired product functionality, and risk of injury, illness or death due to contamination.

- Histoacryl® is made from n-butyl-2-cyanoacrylate (NBCA) and D&C violet No. 2 dye in the coloured version, which are non-metallic and non-conducting materials, therefore it is considered MR safe.

Precautions

- Only product that is liquid in consistency must be used. The condition of the adhesive can be judged with the plastic ampoule unopened.

- Verify that the injection needle is clean and compatible with n-butyl-2-cyanoacrylate and does not trigger polymerisation or degrades with contact. The

injection needle may be degraded if its compatibility has not been confirmed; only those injection needles whose compatibility with cyanoacrylates has been confirmed should be used. Refer to instructions for use regarding any injection needle to be used with Histoacryl® to determine compatibility and prevent injection needle damage

- Care must be taken to avoid early polymerisation within the instruments throughout the application process.

- If an early polymerisation of the Histoacryl® / Ethiodized oil mixture occurs before injection the mixture should be discarded.

- The Histoacryl® /Ethiodized oil mixture shall not be injected if resistance is met during injection. Do not attempt to clear or overcome the resistance by applying increased pressure. If this occurs, determine the cause, and replace the injection needle if necessary.
- Is not recommended to use saline solution or any ionic contrast medium during the procedure. The polymerisation of Histoacryl® starts immediately when it is in contact with an anionic solution.

- Care must be taken to avoid direct contact of instruments, cloths or gloves with the adhesive applied, as they might adhere to the tissue. Instruments stained with Histoacryl® can be cleaned with dimethyl formamide or acetone.

Side effects

The following adverse events have been reported for Histoacryl®: pain, embolisation, stroke, post embolisation syndrome (fever, nausea, and abdominal pain without evidence for infection), fever, end-organ infarction, portal vein thrombosis, splenic vein thrombosis, visceral fistula, aspiration pneumonia, dyspnoea, bacterial peritonitis, bacteremia, bleeding, early or late re-bleeding (that may result in death), liver failure, dysphagia, oesophageal stenosis, cardiac obstruction, intestinal obstruction, bowel infarction, gastric perforation, retrogastric abscess and ulcer.

3. MESH FIXATION

Mode of action

When applied dropwise on the edges of the placed and adapted mesh, the adhesive flows through the mesh

and hardens upon contact with the underlying tissue. The mesh is fixed by the hardened adhesive.

Mode of application

a) Pull open the aluminium pouch and remove the ampoule containing the adhesive in sterile condition. The ampoule containing the adhesive should only be removed from the aluminium pouch immediately prior to use.

b) Hold the ampoule with the cannula pointing upwards and flick it sharply to ensure no adhesive remains trapped within the cannula. Open the ampoule by twisting off the ribbed tip of the cannula keeping it until the end of the procedure. Hold the yellow part of the ampoule between two fingers, keeping it vertical with the thin end pointing upwards while opening the ampoule. This will prevent any loss of Histoacryl® by spillage from the ampoule when opening.

c) Place the mesh in the correct area.

d) Apply some/few small drops of Histoacryl® to the desired fixing sites. The fewer the drops the better the integration will be achieved.

e) For meshes with large pores (more than 1 mm) place the drop on top of the mesh. For meshes with small pores (less than 1 mm) place the drop between the mesh and the tissue.

f) The mesh is held in position for up to 30 seconds to ensure hardening of the glue.

g) Discard any remaining Histoacryl® without contaminating the surroundings. Prior to discarding the ampoule, put back the ribbed tip that was removed when opening the ampoule.

h) Histoacryl® can be used for mesh fixation for the treatment of hernia repair by laparoscopic and open surgery.

Warnings

- Cyanoacrylates, like other fixation materials, may induce the formation of adhesions. Great care must be taken when applying Histoacryl® within the mesh.

- Do not reuse. Reuse might cause an infection hazard for patients and/or users, impaired product functionality, and risk of injury, illness, or death due to contamination.

- Histoacryl® is made from n-butyl-2-cyanoacrylate

(NBCA) and D&C violet No. 2 dye in the coloured version, which are non-metallic and non-conducting materials, therefore it is considered MR safe.

Precautions

- Only product that is liquid in consistency must be used. The condition of the adhesive can be judged with the plastic ampoule unopened.
- Care must be taken to avoid direct contact of instruments, cloths or gloves with the adhesive applied, as they might adhere to the tissue. Instruments stained with Histoacryl® can be cleaned with dimethyl formamide or acetone.
- As general practice, contact of the cannula with the mesh or the tissue should be avoided, in order to minimise the risk of accidental bonding.
- Do not use the yellow end of the ampoule for pushing the mesh after Histoacryl® application. Otherwise, the application end may be adhered to the mesh or the tissue.

Side effects

The following adverse events have been reported for Histoacryl®: pruritus, ecchymosis, hematoma, seroma, surgical site infection, orchitis, hernia recurrence, pain, haemorrhage, mesh infection, bulge, feeling of foreign object, pulmonary migration of glue.

Sterilization

The sterility of Histoacryl® tissue adhesive is ensured by membrane filtration and aseptic filling. Do not re-sterilise.

Histoacryl® is intended for single use. Opened unused ampoules shall be discarded. Opened unused or damaged packs shall be discarded.

Only use Histoacryl® if the package is undamaged. The absence of channels along the sealing of the sterile barrier system must be visually checked as well as the absence of perforations on it. In case of identifying such defects, the device shall be discarded in the recommended manner.

Storage

Histoacryl® shall be stored in its original sealed aluminum pouch in temperatures at or below +22 °C away from moisture and direct light sources.

Histoacryl® shall not be employed after the expiry date

indicated on the label.

Device disposal

Once the treatment is completed the different components of the Histoacryl® must be discarded in special intended containers. It is the responsibility of the user to determine if disposal material is hazardous according to federal, state and local regulations. Dispose of contents and container to comply with applicable local, state, national and international regulation and hospital procedures.

Information to the User/Patient

- In case of a serious incident occurred in relation to Histoacryl®, it should be reported to the Manufacturer and the Competent Authority of the Member State in which the user and/or patient is established.

- When used for skin closure, the polymerised film of Histoacryl® is spontaneously completely detached 4-8 days after application.




- When used for sclerotherapy or mesh fixation, Histoacryl® is intended to degrade slowly in the body. Partial absorption is observed from 30 days being completed at 24 months. Therefore, the expected lifetime of Histoacryl® in the human body is 2 years.


- The clinical staff shall complete the blank spaces of the supplied Implant Card (IC). The meaning of the symbols included in the IC is described in the Section "*Symbol to be used on the label*" of this Instructions for use. The IC shall be given to the patient after the surgical procedure.

- Histoacryl® is a sterile, liquid tissue adhesive made of n-butyl-2-cyanoacrylate (100%). A coloured variant is available which contains n-butyl-2-cyanoacrylate (> 99,985 % w/w) and D&C Violet n°2 0.015 ± 0.001 % w/w)

- Link to the summary of safety and clinical performance on EUDAMED (Art32 MDR). [<https://ec.europa.eu/tools/eudamed/#/screen/home>]

Symbols to be used on the label

	Medical Device
	Dimension(*) (*Ampoule content)
	Catalogue number
	Batch code
	Use-by date
	Date of manufacture
	Sterile unless package is opened or damaged. Method of sterilization: Aseptic Filling
	Single Sterile Barrier System
	Consult Instructions for Use
	Caution
	Upper limit of temperature
	Keep dry
	Keep away from sunlight
	Do not reuse
	Do not use if package is damaged and consult instructions for use
	Unique Device Identifier
	n-butyl-2-cyanoacrylate
	Unfold for more languages

	Patient identification
	Date of implantation
	Health care center or doctor
	Patient information website

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