



Blood Glucose  
Measurement

# Omnitest<sup>®</sup> plus

ACCURACY TEST  
INTERNATIONAL STANDARD ISO 15197:2013

# Omnitest<sup>®</sup> plus

## System Accuracy Evaluation

In accordance with the ISO 15197:2013 standard, new criteria for blood glucose meters have been published.

The aim of this document is to present how Omnitest<sup>®</sup> plus meets these new requirements.

### TEST INFORMATION

#### SYSTEM ACCURACY

The accuracy of Omnitest<sup>®</sup> plus blood glucose monitoring system was assessed by comparing patients' blood glucose results obtained with Omnitest<sup>®</sup> plus with those of a standard laboratory instrument, the YSI 2300 auto analyzer.

1) Test date	November 26 <sup>th</sup> , 2012 - December 12 <sup>th</sup> , 2012
2) Test meter serial number	GKRBMC00034 - GKRBMC00039
3) Test strip lot numbers	B4MK07 (Lot #1), B4MK08 (Lot #2), B4MK09 (Lot #3)
4) Sample numbers	<ul style="list-style-type: none"> <li>▪ 6 x Omnitest<sup>®</sup> plus meter</li> <li>▪ 600 x Omnitest<sup>®</sup> plus test strips (three lots)</li> <li>▪ 1 x YSI 2300 auto analyzer</li> </ul>
5) Standard/guidance documents referenced	<ul style="list-style-type: none"> <li>▪ ISO/DIS 15197:2010 In vitro diagnostic test systems Requirements for blood glucose monitoring systems for self-testing in managing diabetes mellitus (during clinical trial)</li> <li>▪ ISO 15197:2013 In vitro diagnostic test systems Requirements for blood glucose monitoring systems for self-testing in managing diabetes mellitus (for data analysis)</li> <li>▪ CLSI EP09-A2:2004 Method comparison and bias estimating using patient samples</li> </ul>

### SAMPLE DISTRIBUTION

#### SYSTEM ACCURACY

Glucose concentration mmol/L	Percentage of sample	Sample numbers	Preparation of sample
≤ 2.77	5%	5	glycolyzed
> 2.77 - 4.44	15%	15	unaltered
> 4.44 - 6.66	20%	20	unaltered
> 6.66 - 11.10	30%	30	unaltered
> 11.10 - 16.65	15%	15	unaltered
> 16.65 - 22.20	10%	10	unaltered
> 22.20	5%	5	supplemented with glucose
<b>Total</b>	<b>100%</b>	<b>100</b>	

*Distribution of glucose concentrations in samples for system accuracy evaluation*

# System Accuracy Evaluation Difference Plot

## ACCEPTANCE CRITERIA


ACCURACY PLOT – PART 1

Blood glucose concentration & requirements	Tolerance range
< 5.55 mmol/L	± 0.83 mmol/L
≥ 5.55 mmol/L	± 15 %

**Minimum requirement: 95 % of all results within tolerance range**

## TEST RESULTS

SYSTEM ACCURACY PLOT



## DATA ANALYSIS

BIAS DISTRIBUTION ANALYSIS  
COMPARED WITH YSI 2300

### System accuracy results for glucose concentration < 5.55 mmol/L

Strip lot	Within ±0.28 mmol/L	Within ±0.56 mmol/L	Within ±0.83 mmol/L	Combined lots within ± 0.83 mmol/L  <b>99.5 % [185/186]</b>
Lot 1	64.5 % [40/62]	88.7 % [55/62]	100.0 % [62/62]	
Lot 2	72.6 % [45/62]	93.5 % [58/62]	100.0 % [62/62]	
Lot 3	61.3 % [38/62]	91.9 % [57/62]	98.4 % [61/62]	
Combined	66.1 % [123/186]	91.4 % [170/186]	99.5 % [185/186]	

### System accuracy results for glucose concentration ≥ 5.55 mmol/L

Strip lot	Within ±5 %	Within ±10 %	Within ±15 %	Combined lots within ± 15 %  <b>96.4 % [399/414]</b>
Lot 1	55.8 % [77/138]	89.1 % [123/138]	96.4 % [133/138]	
Lot 2	56.5 % [78/138]	89.9 % [124/138]	97.1 % [134/138]	
Lot 3	52.2 % [72/138]	87.0 % [120/138]	95.7 % [132/138]	
Combined	54.8 % [227/414]	88.6 % [367/414]	96.4 % [399/414]	

System accuracy results for glucose concentrations between  
1.7 mmol/L and 31.6 mmol/L  
**Within ± 0.83 mmol/L or ± 15 %**

**97.3 %**  
**[584/600]**

RESULT

97.3 % of results are within ± 0.83 mmol/L or ± 15 % for combined lots.

CONCLUSION

Omnitest® plus exceeds the minimum requirement of 95 % system accuracy.

# System Accuracy Evaluation Consensus Error Grid

## ACCEPTANCE CRITERIA


ACCURACY PLOT – PART 2

Test criteria	Tolerance range
Consensus error grid	99 % of results within zones A and B of the consensus error grid for type 1 diabetes

Zone	Classification
A	No effect on clinical action.
B	Altered clinical action – little or no effect on clinical outcome.
C	Altered clinical action – likely to affect clinical outcome.
D	Altered clinical action – could have significant medical risk.
E	Altered clinical action – could have dangerous consequences.

## TEST RESULTS

CONSENSUS ERROR GRID



## DATA ANALYSIS

CONSENSUS ERROR GRID

Strip lot	Zone A	Zone B	Zone C	Zone D	Zone E	Total
Lot 1	99.5 % [199/200]	0.5 % [1/200]	0 % [0/200]	0 % [0/200]	0 % [0/200]	100 % [200/200]
Lot 2	100 % [200/200]	0.0 % [0/200]	0 % [0/200]	0 % [0/200]	0 % [0/200]	100 % [200/200]
Lot 3	99.5 % [199/200]	0.5 % [1/200]	0 % [0/200]	0 % [0/200]	0 % [0/200]	100 % [200/200]
Combined	99.7 % [598/600]	0.3 % [2/600]	0 % [0/600]	0 % [0/600]	0 % [0/600]	100 % [600/600]

RESULT

CONCLUSION

100 % of results are within zones A and B of the consensus error grid.  
Omnitest® plus exceeds the system accuracy requirement that 99 % of results are within zones A and B.