**Instructions for use – Please read carefully!**

H.E.L.P. NaCl 0.9 %

[1] Product description

Sterile priming and rinsing solution for extracorporeal treatment systems. Only intended for extracorporeal use within the context of H.E.L.P. apheresis.

Steam sterilised and endotoxin-free.

Composition:

1000 ml contains:
Sodium chloride 9.0 g
Water for injection

Which corresponds to:
Na⁺ 154 mmol/l
Cl⁻ 154 mmol/l
Theoretical osmolarity 308 mmol/l

[2] Intended use

The sterile H.E.L.P. NaCl 0.9 % solution is a priming and rinsing solution that is used for the extracorporeal circuit within the context of H.E.L.P. therapy with the Plasmate[®] Futura or Plasmate[®] Secura (summarized below as „Plasmate[®]“).

[3] Indications, contraindications, and side effects

Pay attention to the indications, contraindications and side effects described in the instructions for use of the H.E.L.P. Futura set and in the operating manual of the Plasmate[®].

[4] Warnings and precautions

- Not intended for intravenous infusion!
- The H.E.L.P. NaCl 0.9 % solution must not be used after the expiry date stated on the container and outer packaging
- Only use if the solution is clear and colourless and if the container and connections are undamaged
- Only remove the outer packaging immediately before use.
- Store the solution out of children's reach
- Animal experiments have shown DEHP to be potentially toxic to reproduction. Given the current state of scientific knowledge, a risk (in case of long-term use) especially for premature male babies cannot be completely excluded. As a precaution, the application of medical products containing DEHP should be restricted to short-term use for pregnant women, breastfeeding mothers, infants and children
- FOR SINGLE USE ONLY! DO NOT REUSE!

Re-use of single-use devices creates a potential risk of the patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient.

- DO NOT add any additives
- Discard unused portion

[5] Instructions for use

Remove the outer packaging of the bag and connect it to the H.E.L.P. Futura Set or the tubing of the Plasmate[®] Secura in accordance with the operating manual of the Plasmate[®].

In case of replacement of individual filters within the context of the H.E.L.P. therapy, use H.E.L.P. NaCl 0.9 % solution to flush the new filter(s). The filter replacement procedure is described in the operating manual of the Plasmate[®].

[6] Method of administration and duration of the treatment

Comply with the specifications in the instructions for use of the H.E.L.P. Futura Treatment Set and in the operating manual of the Plasmate[®].

DIN: 02373815

Do not use if packaging is damaged

Store at

B | BRAUN

MANUFACTURER
B. Braun Avitum AG
34209 Melsungen, Germany

CANADIAN DISTRIBUTOR
Chief Medical Supplies Ltd.
411-19th Street S.E.
Calgary, Alberta T2E 6J7

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FABRICANT

B. Braun Avitum AG
34209 Melsungen, Germany

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Chief Medical Supplies Ltd.
411-19th Street S.E.
Calgary, Alberta T2E 6J7

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Steam sterilized and endotoxin-free.

Composition:

1000 mL contains:	
Sodium chloride	9.0 g
Water for injection	

Which corresponds to:

Na ⁺	154 mmol/L
Cl ⁻	154 mmol/L
Theoretical osmolarity	308 mmol/L

[2] Intended use

The sterile H.E.L.P. NaCl 0.9 % solution is a priming and rinsing solution that is used for the extracorporeal circuit within the context of H.E.L.P. therapy with the Plasmate[®] Futura.

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Pay attention to the indications, contraindications, and side effects described in the instructions for use of the H.E.L.P. Apheresis Set and in the operating manual of the Plasmate[®] Futura Apheresis System.

[4] Warnings and precautions

- Not intended for intravenous infusion!
- The H.E.L.P. NaCl 0.9% solution must not be used after the expiry date stated on the container and outer packaging.
- Only use if the solution is clear and colorless and if the container and connections are undamaged.
- Only remove the outer packaging immediately before use.
- Any remaining solution must be discarded.
- Store the solution out of children's reach.

[5] Instructions for use

Remove the outer packaging of the bag and connect it to the H.E.L.P. Futura Set in accordance with the operating manual of the Plasmate[®] Futura Apheresis System.

In case of replacement of individual filters within the context of the H.E.L.P. therapy, use H.E.L.P. NaCl 0.9 %

solution to flush the new filter(s). The filter replacement procedure is described in the operating manual of the Plasmate[®] Futura Apheresis System.

[6] Method of administration and duration of treatment

Follow the specifications in the instructions for use of the H.E.L.P. Apheresis Set and in the operating manual of the Plasmate[®] Futura Apheresis System.

Rx only

Expiry date

Batch number

Manufacture date

Article number

See Instructions for Use

Steam sterilized

Do not re-use

Store at

Do not use if packaging is damaged

Primary packaging:
the container made of PVC
contains DEHP.

B | BRAUN

B. Braun Avitum AG
34209 Melsungen
Germany

US Distributor
B. Braun Medical Inc.
Bethlehem, PA 18018-3524
Made in Germany