

# HYGIENE AUDITING TOOL

## ONBOARDING MOBILE APP

Dr. Thomas-Joerg Hennig  
18 March 2021

# Hygiene Auditing Tool - Global B. Braun APP Drives Digitalization



**Observation Form**

Facility: \_\_\_\_\_ Period Number\*: \_\_\_\_\_ Session Number: \_\_\_\_\_  
 Service: \_\_\_\_\_ Date: (dd/mm/yy) / / Observer: \_\_\_\_\_  
 Ward: \_\_\_\_\_ Start/End time: : / : Page N°: \_\_\_\_\_  
 Department: \_\_\_\_\_ Session duration: (mm) City\*\*: \_\_\_\_\_  
 Country\*\*:



Prof. cat N°	Code	Prof. cat N°	Code	Prof. cat N°	Code	Prof. cat N°	Code
1	<input type="checkbox"/> bef-pat <input type="checkbox"/> bef-asept <input type="checkbox"/> aft-b.f <input type="checkbox"/> aft-pat <input type="checkbox"/> aft.p.surr.	<input type="checkbox"/> HR <input type="checkbox"/> HW <input type="checkbox"/> missed <input type="checkbox"/> gloves	1	<input type="checkbox"/> bef-pat <input type="checkbox"/> bef-asept <input type="checkbox"/> aft-b.f <input type="checkbox"/> aft-pat <input type="checkbox"/> aft.p.surr.	<input type="checkbox"/> HR <input type="checkbox"/> HW <input type="checkbox"/> missed <input type="checkbox"/> gloves	1	<input type="checkbox"/> bef-pat <input type="checkbox"/> bef-asept <input type="checkbox"/> aft-b.f <input type="checkbox"/> aft-pat <input type="checkbox"/> aft.p.surr.



	A	B	C
1			
2			
3			
4			
5			
6			
7			




Time consuming  
Error prone transfer  
Delayed reporting



Time saving  
Results are available in real-time.

\* To be completed by the data manager.  
 \*\* Optional, to be used if appropriate, according to the local needs and regulations.  
 All reasonable precautions have been taken by the World Health Organization to verify the information contained in this document. However, the published material is being distributed without warranty of any kind, either expressed or implied. The responsibility for the interpretation and use of the material lies with the reader. It is to be noted that the World Health Organization has no liability for damages arising from the use of the information. The World Health Organization (WHO), in particular the members of the Infection Control Programme, for their active participation in developing this material. Revised August 2009

# Hygiene Auditing Tool




is a Cloud-Based Customizable Quality Management Tool and a Mobile APP for Data-Entry and Data-Analysis

to Monitor Hand Hygiene Adherence To My 5-Moments & Infection Precautions

for Infection Prevention & Control Professionals and Healthcare Quality Managers who are trained in the World Health Organization hand hygiene observation method.


Sax H, Allegranzi B, Chraïti MN, Boyce J, Larson E, Pittet D. The World Health Organization hand hygiene observation method. Am J Infect Control. 2009 Dec;37(10):827-34. doi: 10.1016/j.ajic.2009.07.003. PMID: 20004812.

# User Role *Institutional Manager* Has Access to Web-Admin-Tool and Mobile APP



## Web-Admin Tool


- Create/edit users
- Create/edit locations
- Set timeline & targets for auditors and check fulfillment rate
- Generate and filter reports
- Download reports, raw data & user lists




## Mobile APP

- Record observations
- Generate and filter reports by location, user and health care worker type


# User Role *Auditor* has Access to the APP, only.




Sees targets (if there are any), fulfillment rate & recorded observations.



Can perform Hand Hygiene Compliance Observations




For the user role *Auditor* the selection of other Auditors is blocked! Data entered by other Auditors cannot be seen.



Generate and filter reports by location, dates & health care worker type

# Main Features Web-Admin-Tool



Reports



Users List



Institutions



Manage HCW



Locations




Profile







Logout?

# Location Setting in the Web-Admin-Tool

## Create Your Own Drop-Down-Menus for the APP





Smart Hospital Group



If no location is defined, the APP will not work!

# Target Setting of No. of Opportunities to be Observed





# Create and Sort Healthcare Worker Groups

**B | BRAUN**  
SHARING EXPERTISE

Icons: Bar chart, List, Calendar, Group, Location, Person, Arrow

## HAT Manager Platform

Healthcare Worker

Add new +

- Nurse x
- Doctor x
- HCA (Healthcare Assistant) x
- Floor Cleaner x
- Technician x**
- Cleaning Contractors x
- Clerical Assistant x
- DSW (Dialysis Support Worker) x
- Spine x

NOTE: Drag to sort items

**B | BRAUN**  
SHARING EXPERTISE

Icons: Bar chart, List, Calendar, Group, Location, Person, Arrow

## HAT Manager Platform

Healthcare Worker


Add new +

- Nurse x
- Technician x**
- Doctor x
- HCA (Healthcare Assistant) x
- Floor Cleaner x
- Cleaning Contractors x
- Clerical Assistant x
- DSW (Dialysis Support Worker) x
- Spine x

NOTE: Drag to sort items

# Main Features Mobile APP

## Language Selection



Currently available language packs

English  
Portuguese (BR)  
Spanish (MX)  
Bahasa Indonesia

# Main Features Mobile APP


## Home Screen

Introduction

Generating Charts

User Profile

Training Videos



Monitoring «My 5-Moments»

Generating Reports

Practice Mode for  
Monitoring «My 5-Moments»

Frequently Asked Questions


# Monitoring «My 5-Moments»

← Back to previous screen

Select HCW-Type


## Monitoring «My 5-Moments»

Moment 4 and 5 cannot be recorded together.  
Select Moment 1-5 by tapping on the arrows.



Link an ACTION to the selected Moments

Location selected on previous screen



Audit summary

Enter Notes  
(max. 150 characters)


Additional Precautions

Move ← to other location

Stop × Observation

Submit → Observation to Server

If the APP is used offline the data will be stored on your mobile device. The next time the APP is used online the data will be submitted to the server and deleted from your mobile device.




## Main Features Mobile APP


### Monitoring My-5-Moments for Hand Hygiene


1. Enter the target number of observations to be met during this observation.
2. Select location
3. Select HCW (Health Care Worker) Type
4. Begin Audit
5. Click on the arrows (Moments) and record action, add notes, precaution measures if indicated then press Submit.
6. Record more hand hygiene opportunities
7. After completing the observations provide educational feedback to the observed HCW

# Main Features Mobile APP

## Monitoring Additional Precautions



1. Tap on Action  - Button
2. Record the required information
3. Drag down to get back to the previous screen.
4. Submit the data




## Main Features Mobile APP


### Generating Real-Time Reports & Charts

1. Click on **Reports** or **Charts**
2. Select the period to review
3. **Filter by Auditor**  
User role **Institutional Manager** can filter by Auditor  
User role **Auditor** this field is blocked.
4. **Filter by HCW** (Health Care Worker)
5. **Filter by location**  
Select different location  
Select sub location  
Select sub sub location
6. Switch to **Reports**  
Every data set can be displayed in full detail  
Location is saved automatically and transferred to **Charts** or vice versa to **Reports**.

# Main Features Mobile APP

## Enter Notes




1. Tap on Action  - Button
2. Enter free text, max. 150 characters




# Main Features Mobile APP

## Practice Mode


To try how to monitor hand hygiene compliance without saving the data to the server.



# Data Security, Data Privacy & Transparency



Same content for the Web-Admin-Tool and the mobile APP. At first-time login Data Privacy and Terms & Conditions must be accepted.



# Transparency

Institutional Managers can download all data that is stored for their Institution

ID	ACTIVE	E	USERNAM	DECRYPTED EMAIL	COMPANY	GROUP NAME	FULL NAME	DATE REGISTERED	EXPIRATION DATE	COUNTRY NAME	COUNTRY CODE	OBSERVATION DATE REGISTERED	TARGET OBSERVATION	OBSERVATIONS MADE	BALANCE START DATE	END DATE
964	1			Test1	B. Braun Medical AG (BBMCH)	Institution Manager	Test7	Oct, 27 2019 01:14 PM	Dec 31, 2020			2019-10-30 13:03:06		0	0	0
962	1			Test2	B. Braun Medical AG (BBMCH)	Institution Manager	Test8	Nov, 05 2019 04:52 PM	Dec 31, 2020			2019-12-16 10:28:09		0	0	0
954	1			Test3	B. Braun Medical AG (BBMCH)	Institution Manager	Test9	Aug, 21 2019 08:55 AM	Dec 31, 2020			2020-05-12 07:29:59		0	0	0
953	1			Test4	B. Braun Medical AG (BBMCH)	auditor	Test10	Aug, 16 2019 09:45 PM	Dec 31, 2020					0	0	0
952	1			Test5	B. Braun Medical AG (BBMCH)	auditor	Test11	Aug, 16 2019 09:44 PM	Dec 31, 2020					0	0	0
951	1			Test6	B. Braun Medical AG (BBMCH)	auditor	Test12	Dec, 16 2019 05:26 PM	Dec 31, 2020			2020-07-07 15:23:52		0	0	0

Date & Time	Auditor	Location	Title	HEALTHCARE WORKER					MOMENT							OCCUPATIONAL			
				1	2	3	4	5	Action	Result	Exposure	GLOVES	Risk	GOWN	MASK	Mask Type	Notes		
2020-01-07 15:39:38	Test	Desk \ Aisle	Nurse	1					rub	PASSED	-	YES	-	-	-				
2020-04-07 12:56:12	Test	Desk \ Window	Doctor		3				rub	PASSED	-	-	-	-					
2020-04-07 12:56:18	Test	Desk \ Window	Doctor					5 missed	FAILED	-	-	-	-						
2020-04-07 12:56:25	Test	Desk \ Window	Doctor		4			wash	PASSED	-	-	-	-						

	Moment 1			Moment 2			Moment 3			Moment 4			Moment 5		
	Passed	Total	%	Passed	Total	%	Passed	Total	%	Passed	Total	%	Passed	Total	%
Nurse (ID 253)	4	5	80%	4	4	100%	5	5	100%	3	5	60%	1	1	100%
Doctor (ID 252)	0	0	0%	1	1	100%	1	2	50%	2	2	100%		2	0%
	0	0	0	0	0	0	0	0	0	0	0	0	0	0	0

# Data Security

The application is served via HTTPS to ensure integrity and confidentiality of data between the user's computer and the site.

The login system flags a user account and locks the account after several failed login attempts. It is currently configured to 5 attempts.

After 15 minutes idle time (no activity on the Web-Admin-Tool) the user is logged off automatically.

Email address of the user field is stored encrypted in the database.

Administrators (Employees from B. Braun) can set up and edit accounts for institutions and users. Recorded HH-compliance data or locations is inaccessible for B. Braun employees.

# Interested?

Please contact  
[hat-service@bbraun.com](mailto:hat-service@bbraun.com)